

381™

Baby Crib

For Rooming-in

MATERNITY



Infant bed Baby Crib 381



Crib trolleys can be horizontally stacked to save place. This is an advantage during transportation and storage.

The status of the information, specifications and illustrations in this publication is indicated by the date of publication. Whilst the greatest care has been taken over its compilation, Merivaara Corp. is not responsible for any errors or omissions. Merivaara Corp. reserves the right to make changes to the technology, features, specifications and design of the equipment and models without notice. Please check the Merivaara website for the latest updated specifications and other information.

All trademarks are Merivaara property unless otherwise stated. Patents granted or pending.


©Merivaara Corp. 6023EN 072012, printed in markprint

For Rooming-in

- Light to handle
- Safe, sturdy construction
- Easy to clean
- Space-saving storage

Infant bed Baby Crib 381 is a handy aid for the nursing of newborns. It has been designed and dimensioned to allow the nurse and the new mother to tend towards the baby in ergonomically correct positions. The large castors make the crib easy to move and they are equipped with brakes to hold the crib safely in place. The crib can be positioned across the bed of the mother, thus allowing convenient access to the baby. The bottom of the frame goes under the bed leaving more clear floor space.

- Frame chrome plated
- Crib of transparent acrylic
- Easy to clean
- Mattress base perforated
- A label pocket at the end
- Storage basket underneath (optional accessory)
- Trendelenburg and reverse Trendelenburg adjustments steplessly
- Castors Ø 100 mm, two castors fitted with brakes



- 381 is classified into the class I according to European Council Directive 93/42/EEC Class I (Annex IX)



Merivaara Corp.
Puustellintie 2, FI-15150 LAHTI, FINLAND
Tel. +358 3 3394 611, fax +358 3 3394 6144
merivaara@merivaara.com
www.merivaara.com