

**NEW**  
Optical  
Design


See beyond expanded field with  
**Panoramic Loupes**






# The Visible Difference<sup>®</sup>

Introducing the first advancement in telescopic optics in 100 years.






Classic Expanded Field Design



Patented Panoramic Design: US# US8928975B2  
(European Patent Pending)



Loupes	3.5X Galilean Standard Field	Micro3.5EF
Viewable Area	20 cm <sup>2</sup>	44 cm <sup>2</sup>
Field Diameter	5 cm	7.5 cm
Angle of View	8°	12°
<div>  <b>Viewing Field</b> 18" working distance</div>		

# Panoramic Loupes

See the bigger picture with the new patented panoramic field design, doubling the viewable area over prismatic expanded field optics. Designs for Vision's new panoramic loupes provide precision magnification with an unprecedented field of view.



See the Visible Difference today with a 45 day no obligation trial!



**3.5X Classic Expanded Field**

50 cm<sup>2</sup>

8 cm

14°





**NEW 3.5X Panoramic Field**

100 cm<sup>2</sup>

11.1 x 9 cm

20°





# HDi™

## High Definition Imaging

Patent Pending



without HDi



with HDi

Advanced photonic design provides  
uniform light distribution with  
maximum intensity



## LED DayLite® UltraMini HDi

*Small, lightweight and modular*

Now available with a Panoramic Spot, illuminating  
the entire viewing field of the panoramic loupes.



**WE'VE EXPANDED** – Note our NEW ADDRESS

4000 Veterans Memorial Hwy. Bohemia, NY 11716-1024

1.800.345.4009

1.631.585.3300

info@dvimail.com

DesignsForVision.com

