



Instrument and apparatus trolleys
multi-purpose and transport trolleys

INFORMATION

505
E

Instrument, apparatus and multi-purpose trolleys

comply with the DIN standard for packaging of sterile goods.



Instrument trolleys are available in height 800, 900 and 1000 mm

They consist of a sturdy chrome nickel steel tubular frame and are mobile on 4 antistatic swivel castors, 75 mm diameter, with bumpers. Two of the castors are fitted with brakes. Removable chrome nickel steel shelves are either inserted or placed on top. Optional extras include a plastic guard and a drawer in closed housing.

Apparatus trolleys, 650 mm high

Guarantee an optimum working height for locating equipment. For the technical specifications, see the instrument trolleys.



Instrument and apparatus trolleys may be provided with twin wheel castors to endure MRI compatibility.

BECKER SYSTEMS for:

wards
out-patient departments
operating rooms
treatment rooms
sterilization

BECKER



Multi-purpose trolleys, 865 mm high feature a **push handle** and **swivel castors**, diameter **125 mm** for greater manoeuvrability. In all other respects, the design of these trolleys is identical to that of the “instrument trolleys”.



Transport trolleys, 940 or 1260 mm high, are provided **with shelves firmly welded in between two circular tube push handles**, featuring a high rim with large radii all around to avoid soiling. Mobile on 4 castors with a diameter of 125 mm. 2 of the castors are fitted with brakes. Plastic bumpers are fitted to prevent damage to the trolley.

Material

All trolleys and tables as well as the chrome-nickel steel accessories are mainly made of high-grade chrome nickel steel sheets and tubing, material AISI 304 (1.4301), surface polished grain 240/320.

For specific-purpose use, we also offer **anaesthetic table, sterilization tray trolleys, pallet trolleys, and universal trolleys**; they also comply with the DIN standard for packaging of sterile goods.

Instrument trolleys, 800 mm, 900 mm or 1000 mm high and anaesthetic table

complying with DIN standard for packaging of sterile goods. Made of high-grade chrome nickel steel sheets and tubing, material AISI 304 (1.4301), surface polished grain 240/320.

Mobile on 4 antistatic swivel castors, diameter 75 mm, double-row ball bearings, anti-thread mechanism, bumpers, two of the castors provided with brakes.

		Height: 800 mm		
	DESCRIPTION	DIMENSIONS (WxDxH)	ORDER NO.	
	Instrument trolley 2 removable chrome nickel steel shelves	625 x 425 x 800 mm 720 x 625 x 800 mm 940 x 625 x 800 mm	150001 150003 150006	
	Instrument trolley 2 removable chrome nickel steel shelves, top with plastic guard on three sides	625 x 425 x 800 mm 720 x 625 x 800 mm 940 x 625 x 800 mm Height incl. guard 875 mm	150011 150013 150016	
	Instrument trolley 2 removable chrome nickel steel shelves, drawer with closed housing	625 x 425 x 800 mm 720 x 625 x 800 mm 940 x 625 x 800 mm	150021 150023 150026	
	Instrument trolley 2 removable chrome nickel steel shelves, drawer with closed housing, top with plastic guard on three sides	625 x 425 x 800 mm 720 x 625 x 800 mm 940 x 625 x 800 mm Height incl. guard 875 mm	150031 150033 150036	
	Instrument trolley with chrome nickel steel shelves; top shelf: smooth and rimless lower shelf: removable	625 x 425 x 800 mm 720 x 625 x 800 mm 940 x 625 x 800 mm	150041 150043 150046	
	Anaesthetic table with removable chrome nickel steel basin 500 x 300 x 100 mm and a removable chrome nickel steel shelf, below	625 x 425 x 800 mm	150080	

Height: 900 mm
Height: 1000 mm

Height: 900 mm			Height: 1000 mm		
DIMENSIONS (WxDxH)	ORDER NO.		DIMENSIONS (WxDxH)	ORDER NO.	
625 x 425 x 900 mm 720 x 625 x 900 mm 940 x 625 x 900 mm	150101 150103 150106		625 x 425 x 1000 mm 720 x 625 x 1000 mm 940 x 625 x 1000 mm	150201 150203 150206	
625 x 425 x 900 mm 720 x 625 x 900 mm 940 x 625 x 900 mm Height incl. guard 975 mm	150111 150113 150116		625 x 425 x 1000 mm 720 x 625 x 1000 mm 940 x 625 x 1000 mm Height incl. guard 1075 mm	150211 150213 150216	
625 x 425 x 900 mm 720 x 625 x 900 mm 940 x 625 x 900 mm	150121 150123 150126		625 x 425 x 1000 mm 720 x 625 x 1000 mm 940 x 625 x 1000 mm	150221 150223 150226	
625 x 425 x 900 mm 720 x 625 x 900 mm 940 x 625 x 900 mm Height incl. guard 975 mm	150131 150133 150136		625 x 425 x 1000 mm 720 x 625 x 1000 mm 940 x 625 x 1000 mm Height incl. guard 1075 mm	150231 150233 150236	
625 x 425 x 900 mm 720 x 625 x 900 mm 940 x 625 x 900 mm	150141 150143 150146		625 x 425 x 1000 mm 720 x 625 x 1000 mm 940 x 625 x 1000 mm	150241 150243 150246	

Apparatus trolleys 650 mm high and sterilization tray trolleys

complying with DIN standard for packaging of sterile goods. Made of high-grade chrome nickel steel sheets and tubing, material AISI 304 (1.4301), surface polished grain 240/320.

Mobile on 4 antistatic swivel castors, diameter 75 mm, double-row ball bearings, anti-thread mechanism, bumpers, two of the castors provided with brakes.






	DESCRIPTION	DIMENSIONS (WxDxH)	ORDER NO.	
	Apparatus trolley 2 removable chrome nickel steel shelves	625 x 425 x 650 mm 720 x 625 x 650 mm	150301 150303	
	Apparatus trolley 2 removable chrome nickel steel shelves, top with plastic guard on three sides	625 x 425 x 650 mm 720 x 625 x 650 mm Height incl. guard 725 mm	150311 150313	
	Apparatus trolley 2 removable chrome nickel steel shelves, drawer with closed housing	625 x 425 x 650 mm 720 x 625 x 650 mm	150321 150323	
	Apparatus trolley 2 removable chrome nickel steel shelves, drawer with closed housing, top with plastic guard on three sides	625 x 425 x 650 mm 720 x 625 x 650 mm Height incl. guard 725 mm	150331 150333	
 for wire trays size 500 x 300 mm	Sterilization tray trolley top with chrome nickel steel frame for wire trays and removable chrome nickel steel, lower shelf (supplied without wire trays) Sterilization tray trolley chrome nickel steel frame for 2 wire trays on top and 2 wire trays below (supplied without wire trays) Perforated wire tray with folding handles	720 x 625 x 800 mm 720 x 625 x 800 mm 480 x 250 x 50 mm	150380 150383 120665	

MULTI-PURPOSE TROLLEYS

Multi-purpose trolleys 865 mm high

complying with DIN standard for packaging of sterile goods. Made of high-grade chrome nickel steel sheets and tubing, material AISI 304 (1.4301), surface polished grain 240/320.



For greater manoeuvrability with a chrome nickel steel push handle and 4 antistatic swivel castors with a diameter of 125 mm, double-row ball bearings, anti-thread mechanism, bumpers, two of the castors provided with brakes.

	DESCRIPTION	DIMENSIONS (WxDxH)	ORDER NO.	
	<p>Multi-purpose trolley 2 removable chrome nickel steel shelves and push handle</p> <p>Shelf size: 720 x 625 mm Shelf size: 940 x 625 mm</p>	<p>830 x 625 x 865 mm 1050 x 625 x 865 mm</p>	<p>150403 150406</p>	
	<p>Multi-purpose trolley 2 removable chrome nickel steel shelves, push handle and top with plastic guard on three sides</p> <p>Shelf size: 720 x 625 mm Shelf size: 940 x 625 mm</p>	<p>Height incl. guard 940 mm</p> <p>830 x 625 x 865 mm 1050 x 625 x 865 mm</p>	<p>150413 150416</p>	
	<p>Multi-purpose trolley 2 removable chrome nickel steel shelves, push handle and drawer with closed housing</p> <p>Shelf size: 720 x 625 mm Shelf size: 940 x 625 mm</p>	<p>830 x 625 x 865 mm 1050 x 625 x 865 mm</p>	<p>150423 150426</p>	
	<p>Multi-purpose trolley 2 removable chrome nickel steel shelves, push handle, drawer with closed housing and top with plastic guard on three sides</p> <p>Shelf size: 720 x 625 mm Shelf size: 940 x 625 mm</p>	<p>Height incl. guard 940 mm</p> <p>830 x 625 x 865 mm 1050 x 625 x 865 mm</p>	<p>150433 150436</p>	
	<p>Multi-purpose trolley with chrome nickel steel shelves; top shelf: smooth and rimless lower shelf: removable push handle</p> <p>Shelf size: 720 x 625 mm Shelf size: 940 x 625 mm</p>	<p>830 x 625 x 865 mm 1050 x 625 x 865 mm</p>	<p>150443 150446</p>	
	<p>Optional extra: 3rd removable chrome nickel steel shelf (welded in between upper and lower shelf) (Only possible when ordered with trolley.)</p>	<p>720 x 625 mm 940 x 625 mm</p>	<p>150453 150456</p>	

Pallet trolleys and universal trolleys

complying with DIN standard for packaging of sterile goods. Made of high-grade chrome nickel steel sheets and tubing, material AISI 304 (1.4301), surface polished grain 240/320.




Mobile on 4 antistatic swivel castors, diameter 125 mm, double-row ball bearings, anti-thread mechanism, bumpers, two of the castors provided with brakes.

	DESCRIPTION	DIMENSIONS (WxDxH)	ORDER NO.	
	<p>Pallet trolley with removable chrome nickel steel shelf and handle</p> <p>Shelf size: 720 x 625 mm Shelf size: 940 x 625 mm</p> <p>Height top edge shelf level: 290 mm</p>	<p>830 x 625 x 860 mm 1050 x 625 x 860 mm</p>	<p>150503 150506</p>	
	<p>Universal trolley 3 removable chrome nickel steel shelves and push handle</p> <p>Shelf size: 720 x 625 mm Shelf size: 940 x 625 mm</p>	<p>830 x 625 x 1050 mm 1050 x 625 x 1050 mm</p>	<p>150513 150516</p>	

Transport trolleys







Tubular frame construction of chrome nickel steel, diameter 25 mm, with 2 integrated handles and welded, rimmed chrome nickel steel shelves.

Mobile on 4 antistatic swivel castors, diameter 125 mm, double-row ball bearings, anti-thread mechanism, bumpers, two of the castors provided with brakes.

	<p>Transport trolley 2 chrome nickel steel shelves</p> <p>Shelf size: 800 x 500 mm Shelf size: 1000 x 600 mm</p>	<p>900 x 600 x 940 mm 1100 x 700 x 940 mm</p>	<p>150550 150555</p>	
	<p>Transport trolley 3 chrome nickel steel shelves</p> <p>Shelf size: 800 x 500 mm Shelf size: 1000 x 600 mm</p>	<p>900 x 600 x 940 mm 1100 x 700 x 940 mm</p>	<p>150560 150565</p>	
	<p>Transport trolley 4 chrome nickel steel shelves</p> <p>Shelf size: 800 x 500 mm Shelf size: 1000 x 600 mm</p>	<p>900 x 600 x 1260 mm 1100 x 700 x 1260 mm</p>	<p>150570 150575</p>	

ACCESSORIES

	DESCRIPTION	DIMENSIONS (WxDxH)	ORDER NO.	
	Dispenser for syringes 3 tilting bins, chrome nickel steel tablet and fixing bracket	600 x 200 x 360 mm	150700	
	Dispenser for syringes 3 bins First-in / First-out-system, chrome nickel steel tablet and fixing bracket	600 x 200 x 360 mm	150705	
	Dispenser for syringes and cannulas 3 big and 5 small tilting bins with chrome nickel steel fixing bracket	600 x 200 x 520 mm	150710	
	Dispenser for syringes and cannulas 3 big and 5 small bins First-in / First-out-system, with chrome nickel steel fixing bracket	600 x 200 x 520 mm	150715	
	Dispenser set for syringes and cannulas 3 big and 5 small tilting bins with combiset and chrome nickel steel fixing bracket	750 x 200 x 520 mm	150720	
	Dispenser set for syringes and cannulas 3 big and 5 small bins First-in / First-out-system, with combiset and chrome nickel steel fixing bracket	750 x 200 x 520 mm	150725	
 	Holder with bowl chrome nickel steel	352 x 325 x 100 mm	150730	
	Lid for bowl 352 x 325 mm chrome nickel steel		110268	
	Holder with bucket 10 l chrome nickel steel	Ø 300 x 250 mm	150740	
	Lid for bucket 10 l chrome nickel steel	Ø 300 mm	110273	
	Ampoule opener with collecting bin and diamond file plus holder		150750	
	Polypropylene disposal containers set of 10 with chrome nickel steel holder	Ø 140 x 145 mm	150760	

	DESCRIPTION	DIMENSIONS (WxDxH)	ORDER NO.	
	<p>Plastic waste container with lid and 1 roll of 48 plastic bags (15 l)</p>	300 x 210 x 440 mm	150800	
	<p>Chrome nickel steel waste container with lid and 1 roll of 48 plastic bags (15 l)</p>	285 x 180 x 400 mm	150805	
	<p>Spare roll of 48 plastic bags (15 l)</p>		110228	
	<p>Waste disposal set A with plastic waste container</p>		150810	
	<p>Waste disposal set B with chrome nickel steel waste container</p>		150815	
	<p>Chrome nickel steel holder with polypropylene container and closed lid and slotted lid</p>	300 x 225 x 130 mm	150820 150825	
	<p>Polypropylene (PP) container sterilizable up to 134°C, for used instruments, syringes, cannulas and catheters with closed lid with slotted lid</p>	300 x 225 x 130 mm	120680 120685	
	<p>Wire sieve tray chrome nickel steel, closely meshed with carrying handle, four feet, drip and stacking frame, for polypropylene containers</p>	245 x 175 x 100 mm	120690	

ACCESSORIES

	DESCRIPTION	DIMENSIONS (WxDxH)	ORDER NO.	
	Holder with swing-out tray chrome nickel steel	235 x 190 x 40 mm	150830	
	Chrome nickel steel tray	235 x 190 x 40 mm	110288	
	Catheter receptacle 3 compartments, chrome nickel steel, to hang on		150840	
	Infusion bottle holder with 4 hooks chrome nickel steel (Adjustable length 700 mm)		150850	
	5-socket connector strip with power indicator light and 2.2 m cable		150860	
	Low-voltage halogen spot lamp 3 m cable, 230 V ac with mounting clamp		150870	
	Surcharge for castor-set (set of four) plastic "thin wheel castors", dia. 75 mm, 2 castors with brakes, bumpers, MRI compatibility (Suitable for all trolleys with castors of 75 mm diameter)		150880	

BECKER



Becker quality steel furniture and equipment have been in use in hospitals and medical centers for decades. Optimum use of available space is achieved through an extensive line of programmes and sizes.

Special designs and individual choice of colors as well as custom-built systems offer an enormous variety of options to the customer.

Metallwarenfabrik
Walter H. Becker GmbH
Anzenkirchener Straße 4
D-84371 Triftern

Tel.: +49 (0) 85 62 / 96 05-0
Fax: +49 (0) 85 62 / 24 34
E-mail: becker.triftern@t-online.de
Internet: www.becker-triftern.de

BECKER SYSTEMS for:

wards
out-patient departments
operating rooms
treatment rooms
sterilization