



**BECKER**

Movable Equipment  
for OR, Outpatient and  
Sterilisation Departments

**INFORMATION**

**310  
E**




# Movable Equipment for OR, Outpatient and Sterilisation Departments

## Movable Equipment for OR, Outpatient and Sterilisation Departments

made of quality chrome nickel steel (CNS) tubes and plates, alloy AISI 304 (1.4301), with polished surface to grain 240/320. All castors and legs conduct electricity.

ILLUSTRATION	DESCRIPTION	ART. NO.
	<p><b>Mayo instrument table</b> height manually adjustable 850 – 1280 mm, 4 swivel castors 75 mm Ø, 2 with brakes</p> <p>Tray size: 600 x 400 mm Tray size: 700 x 500 mm</p>	<p>130010 130013</p>
	<p><b>Oil-pump operated Mayo instrument table</b>, deep-drawn tray with rounded corners, can be swiveled around 360°, lockable in all positions, height adjustable oil hydraulically 925 – 1310 mm, 3 swivel castors 50 mm Ø</p> <p>Tray size: 600 x 400 mm Tray size: 700 x 500 mm Tray size: 850 x 600 mm</p> <p>height adjustable oil hydraulically 950 – 1335 mm, 3 swivel castors 75 mm Ø, 2 with brakes</p> <p>Tray size: 600 x 400 mm Tray size: 700 x 500 mm Tray size: 850 x 600 mm</p>	<p>130020 130023 130026  130021 130024 130027</p>
	<p><b>Oil-pump operated Mayo instrument table</b> deep-drawn tray with rounded corners, top tray can be swiveled 360° and locked in all positions, tray format 1220 x 620 mm height adjustable oil hydraulically 985 – 1370 mm, 4 swivel castors 100 mm Ø, 2 with brakes</p> <p>with 1 tray (top) with 2 trays (top and bottom)</p>	<p>130031 130036</p>
	<p><b>Arm rest</b> height manually adjustable 790 bis 1220 mm, star base 600 mm Ø, 5 swivel castors 75 mm Ø, 2 with brakes</p> <p>with fixed CNS table top 450 x 325 mm (without cushion) with adjustable CNS table top 450 x 325 mm (without cushion)</p> <p><b>Cushion</b> 450 x 325 mm, 40 mm high, black, electrically conductive, for arm rest no. 130050 or 130055</p>	<p>130050 130055 130058</p>
	<p><b>Hand and arm rest</b> height manually adjustable 680 – 980 mm, table top 470 x 450 mm radiolucent, ball-type swivel joint tilts in all directions (without cushion), 4 swivel castors 75 mm Ø, 2 with brakes</p> <p><b>Cushion</b> 470 x 450 mm, 40 mm high, black, electrically conductive, radiolucent, for hand and arm rest no. 130060</p>	<p>130060  130068</p>

ILLUSTRATION	DESCRIPTION	ART. NO.
	<p><b>Kick bucket without push handle</b>, height 425 mm, 4 swivel castors 50 mm Ø, star base 450 mm Ø</p> <p>with chrome nickel bucket 15 l with plastic bucket, white, 10 l</p>	<p>130110 130115</p>
	<p><b>Kick bucket with push handle</b>, height 825 mm, 4 swivel castors 50 mm Ø, star base 450 mm Ø</p> <p>with chrome nickel bucket 15 l with plastic bucket, white, 10 l</p> <p><b>Chrome nickel steel lid</b> for chrome nickel steel bucket 15 l <b>Chrome nickel steel bucket</b> 15 l <b>Plastic bucket</b>, white, 10 l</p>	<p>130120 130125 130127 130128 130129</p>
	<p><b>Single bowl stand</b> 1 chrome nickel steel bowl 8 l, push handle</p> <p>with 4 swivel castors 50 mm Ø [500 x 500 mm, height 800/965 mm] with 5 swivel castors 75 mm Ø, 2 with brakes [650 mm Ø, height 800/965 mm]</p> <p><b>Spare bowl</b> 8 l, chrome nickel steel [220/360 mm Ø, height 150 mm]</p>	<p>130140 130141 130148</p>
 <p>130146                      130145</p>	<p><b>Double bowl stand</b> 2 chrome nickel steel bowls 8 l, push handle</p> <p>with 4 swivel castors 50 mm Ø [500 x 790 mm, height 800/965 mm] with 5 swivel castors 75 mm Ø, 2 with brakes [650 mm Ø, above 760 mm width, height 800/965 mm]</p> <p><b>Bowl stand with heating</b></p>	<p>130145 130146 130150</p>
 <p>Ø 400</p>	<p><b>Linen hamper</b> with linen bag (ring 400 mm Ø) and push handle</p> <p>with 4 swivel castors 50 mm Ø [550 x 550 mm, height 800/965 mm] with 5 swivel castors 75 mm Ø, 2 with brakes [650 mm Ø, height 800/965 mm]</p> <p><b>Spare linen bag</b></p>	<p>130160 130161 130168</p>
 <p>130170 + 130178</p>	<p><b>Instrument receptacle</b> with push handle (without wire tray)</p> <p>with 4 swivel castors 50 mm Ø [550 x 550 mm, height 830/965 mm] with 5 swivel castors 75 mm Ø, 2 with brakes [650 mm Ø, height 830/965 mm ]</p> <p><b>Perforated wire tray with folding handles</b>, 240 x 250 x 50 mm <b>Perforated wire tray with folding handles</b>, 480 x 250 x 50 mm</p>	<p>130170 130171 130177 130178</p>



# Movable Equipment for OR, Outpatient and Sterilisation Departments














ILLUSTRATION	DESCRIPTION	ART. NO.
	<p><b>Surgeon's single step</b> removable, grooved non-slip plate 600 x 350 mm</p> <p>height 130 mm height 220 mm height 300 mm</p>	<p>130200 130203 130206</p>
	<p><b>Surgeon's double step</b> removable, grooved non-slip plates 600 x 350 mm [600 x 650 mm, height 220/440 mm]</p>	<p>130210</p>
	<p><b>Wall-mounted system to store 4 surgeon's steps</b> pair of rails for wall mounting with hooks to hold max. 4 surgeon's single steps no. 130200, 130203 and/or 130206; (without surgeon's steps) [600 x 430 mm, height 1450 mm]</p>	<p>130220</p>
 <p>130243      130246</p>	<p><b>Infusion stands</b> height adjustable 1370 to 2450 mm, 4 hooks, 5 swivel castors 75 mm Ø, 2 with brakes, star base 650 mm Ø</p> <p>with PP drip catcher "Model Ulm" without drip catcher, with push handle</p>	<p>130240 130246 130243</p>
	<p><b>Infusion stand</b> with one-hand height adjustment from 1350 to 2150 mm, with drip catcher, 4 hooks, 5 castors 75 mm Ø, 2 with brakes, star base 635 mm Ø</p> <p><b>Infusion stand</b> for fastening to a standard equipment rail (See Register/Brochure 1)</p>	<p>130250</p>

ILLUSTRATION	DESCRIPTION	ART. NO.
	<p><b>Linen and waste collector</b> with push handle, insertable bottom plate, 4 swivel castors 75 mm Ø, 2 with brakes (without lid and bags)</p> <p>Single [465 x 490 mm, height 900 mm] Double [770 x 490 mm, height 900 mm]</p>	<p>130300 130305</p>
	<p><b>Linen and waste trolley</b> with Becker clips, removable bottom plate, 4 swivel castors 75 mm Ø, (without bags)</p> <p>Single [425 x 425 mm, height 900 mm] Double [970 x 425 mm, height 900 mm]</p> <p><b>Chrome nickel steel lid</b> for linen and waste collectors or trolleys, Art. No.130300, 130305, 130311 and 130314</p> <p><b>Spare clip</b></p>	<p>130311 130314 130308 130317</p>
	<p><b>Hook bar for OR aprons</b> for wall-mounting, double hooks, movable</p> <p>with 3 double hooks, 500 mm long with 4 double hooks, 600 mm long with 5 double hooks, 800 mm long with 6 double hooks, 1000 mm long with 8 double hooks, 1200 mm long</p> <p><b>Spare double hooks</b> for hanging</p>	<p>130360 130362 130364 130366 130368 130378</p>
	<p><b>Operating room shoe rack</b> with 6 hanging wire baskets for 24 pairs of shoes</p> <p>mobile on 4 swivel castors 75 mm Ø, 2 with brakes [910 x 600 mm, height 1600 mm]</p> <p>for wall-mounting with hook bars (2) [910 x 280 mm, height 1400 mm]</p> <p><b>Spare CNS tray</b> [900 x 265 mm, height 35/60 mm]</p>	<p>130410 130415 130418</p>
 <p>130431+5 x 410</p>	<p><b>Trolley for linen &amp; surgical drapes</b> with 3 hanging wire baskets [845 x 560 x 400/180 mm high] (without wire dividers)</p> <p>on 4 swivel castors 75 mm Ø, 2 with brakes [980 x 590 mm, height 1635 mm]</p> <p>on 4 swivel castors 100 mm Ø, 2 with brakes [980 x 590 mm, height 1660 mm]</p> <p><b>Wire divider</b> for baskets 560 mm deep</p>	<p>130430 130431 410</p>



# Movable Equipment for OR, Outpatient and Sterilisation Departments

ILLUSTRATION	DESCRIPTION	ART. NO.
	<p><b>Paper rack</b> for 4 different types of paper up to 1300 mm long, 4 swivel castors 100 mm Ø, 2 with brakes [1400 x 650 mm, height 920 mm]</p>	130450
	<p><b>Paper rack, inclined</b> holds 5 different types of paper up to 1300 mm long 4 swivel castors 100 mm Ø, 2 with brakes [1380 x 580 mm, height 955/1155 mm]</p>	130455
	<p><b>Container transport trolley</b> for sterile material containers 600/300 mm, 4 swivel castors 125 mm Ø, 2 with brakes</p> <p>one-sided with 6 chrome nickel steel trays [650 x 450 mm, height 1540 mm]</p> <p>two-sided with 12 chrome nickel steel trays [650 x 750 mm, height 1540 mm]</p> <p><b>Spare CNS tray</b> [620 x 320 mm, height 175 mm]</p>	<p>130470</p> <p>130475</p> <p>130478</p>
 <p>130525</p>	<p><b>Shelf trolley</b> with fix welded shelves, 4 swivel castors 125 mm Ø, 2 with brakes</p> <p>with 3 non-perforated shelves with 3 perforated shelves [1340 x 670 mm, height 1170 mm, clear distance between shelves 395 mm]</p> <p>with 4 non-perforated shelves with 4 perforated shelves [1340 x 670 mm, height 1600 mm, clear distance between shelves 395 mm]</p> <p>with 5 non-perforated shelves with 5 perforated shelves [1340 x 670 mm, height 1600 mm, clear distance between shelves 288 mm]</p>	<p>130500</p> <p>130505</p> <p>130510</p> <p>130515</p> <p>130520</p> <p>130525</p>
	<p><b>Shelf linen trolley</b> with 4 fixed welded non-perforated shelves, 9 compartment dividers and 6 sets of side dividers, 4 swivel castors 125 mm Ø, 2 with brakes [1340 x 670 mm, height 1600 mm]</p>	130550

# OR Trolleys M9 and M11



Made of powder-coated sheet steel, in a colour selected by the customer from the Becker colour list. Self-contained cabinet construction with double-walled sides and cover plate. 2 welded-on handles. 2 lockable swinging doors opening 270°, doors can be fixated in open position. Sturdy tubular frame construction (with wall protection strips) on 4 swivel castors, 125 mm Ø, 2 with brakes.

Available in two sizes: **M9 trolley**, external dimensions: 1030 x 760 mm, 1450 mm high  
**M11 trolley**, external dimensions: 1230 x 760 mm, 1450 mm high

ILLUSTRATION	DESCRIPTION	ART. NO.
	<p><b>OR modular trolley M9</b> with middle dividing wall to hold inside panels with slots for modules</p> <p>Details on interior fittings supplied on request.</p>	<p>BM09 100</p>
	<p><b>OR modular trolley M9 fitted</b> with two pairs of inside panels with slots for modules 6 module baskets 600 x 400 mm, 200 mm high and 4 module baskets 600 x 400 mm, 100 mm high</p> <p>without dividers – see Brochure 1: Anaesthetic and Dressing Trolleys, page 1.21</p>	<p>BM09 105</p>
 <p style="text-align: right;">BM11 205</p>	<p><b>OR trolley for drapes</b> fitted with one chrome-nickel-steel basket and one interior drawer, 100 mm high with chrome-nickel-steel wire railing</p> <p>Type M9 (width 1030 mm) Type M11 (width 1230 mm)</p> <p><b>Wire divider</b> for basket, chrome nickel steel <b>Wire divider</b> for railing, chrome nickel steel</p>	<p>BM09 205 BM11 205 410 415</p>
 <p style="text-align: right;">BM11 305</p>	<p><b>OR trolley for disposables</b> fitted with 3 interior drawers, 100 mm high, and 2 drawers, 185 mm high, (without dividers)</p> <p>Type M9 (width 1030 mm) Type M11 (width 1230 mm)</p>	<p>BM09 305 BM11 305</p>
 <p style="text-align: right;">BM11 350</p>	<p><b>OR trolley for sterile implants</b> fitted with 7 interior drawers, 100 mm high, (without dividers)</p> <p>Type M9 (width 1030 mm) Type M11 (width 1230 mm)</p>	<p>BM09 350 BM11 350</p>

*For information on drawer dividers, refer to the brochures M9 and M11 as well as the brochures U21 and U10. More information on the trolleys with individual interior fittings can be found in the brochures M9 and M11.*

# BECKER



For decades Becker quality steel furniture and equipment have been a fixed item in medical technology in clinics and hospitals.

Design, individual choice of colours as well as an extensive range of items and sizes allow you to make optimum use of your available space. The resulting concept is customized to your needs and logical through and through.



Metallwarenfabrik  
Walter H. Becker GmbH  
Anzenkirchener Strasse 4  
84371 Triftern, Germany

Phone +49(0)8562 9605-0  
Fax +49(0)8562 9605-64  
info@becker-triftern.de  
www.becker-triftern.de

### Becker equipment:

- for the ward
- for the operating theatre
- for the central sterilisation unit
- for all functional rooms