



BECKER

Swivel Stools,
Mobile Swivel Seats
& Chairs for OR

INFORMATION

4.1
E

2015



Static and Mobile Swivel Stools for OR

- Sturdy five star base (Ø 580 mm) made of chrome nickel steel AISI 304 (1.4301) electropolished
- Concealed precision spindle for infinitely variable height adjustment
- Cushioned seat (black) conductive, Ø 340 mm

ILLUSTRATION	DESCRIPTION	ART. NO.
	<p>Swivel stool, static, with chrome nickel steel seat, height adjustable 435 – 650 mm</p>	140010
	<p>Swivel stool, static, with 40 mm high cushioned seat, height adjustable 465 – 680 mm</p>	140013
	<p>Swivel stool, static, with 65 mm high cushioned seat, height adjustable 490 – 705 mm</p>	140016
	<p>Mobile swivel stool with chrome nickel steel seat, height adjustable 435 – 650 mm 5 conductive safety castors (Ø 50 mm) with internal weight-activated brake option, free rolling under load, appropriate braking action when not under load</p>	140030
	<p>Mobile swivel stool with 40 mm high cushioned seat, height adjustable 465 – 680 mm 5 conductive safety castors (Ø 50 mm) with internal weight-activated brake option, free rolling under load, appropriate braking action when not under load</p>	140033
	<p>Mobile swivel stool with 65 mm high cushioned seat, height adjustable 490 – 705 mm 5 conductive safety castors (Ø 50 mm) with internal weight-activated brake option, free rolling under load, appropriate braking action when not under load</p>	140036
	<p>Back rest for static and mobile swivel stools with cushioned seat (no retrofitting possible)</p>	140050



Mobile Swivel Seats for OR

- Sturdy five-star base (Ø 610 mm) made of powder-coated aluminium
- 5 conductive safety castors (Ø 50 mm) with internal weight-activated brake option, free rolling under load, appropriate braking action when not under load
- Infinitely variable pneumatic height adjustment 530 – 710 mm, operated by footring that is accessible from any seating position
- Cushioned seat (black), conductive, Ø 360 mm, 65 mm high
- Back rest height adjustable for an optimum ergonomic working position

ILLUSTRATION	DESCRIPTION	ART. NO.
	<p>Mobile swivel seat</p> <p>Mobile swivel seat with back rest</p> <p>Star base and column powder-coated colour White Aluminium (RAL 9006)</p>	<p>140100</p> <p>140105</p>
	<p>Mobile swivel seat</p> <p>Mobile swivel seat with back rest</p> <p>Star base and column powder-coated in a single colour selected by the customer from the Becker colour list</p>	<p>140120</p> <p>140125</p>

Mobile Swivel Chairs for OR

- Sturdy non-tilting base plate (Ø 470 mm) with chrome nickel steel surface AISI 304 (1.4301)
- 5 conductive safety castors (Ø 50 mm) with internal weight-activated brake option, free rolling under load, appropriate braking action when not under load
- Chair additional lockable via foot brake
- Foot-operated pneumatic height adjustment 600 – 750 mm, infinitely variable
- Cushioned seat (black), conductive
- Height and tilt angle of back rest adjustable for an optimum ergonomic working position

	<p>Mobile swivel chair</p> <p>with round cushioned seat (Ø 360 mm, 65 mm high)</p> <p>Mobile swivel chair</p> <p>with round cushioned seat (Ø 360 mm, 65 mm high) and back rest</p>	<p>140210</p> <p>140215</p>
	<p>Mobile swivel chair</p> <p>with cushioned saddle seat (375 x 420 mm)</p> <p>Mobile swivel chair</p> <p>with cushioned saddle seat (375 x 420 mm) and back rest</p>	<p>140220</p> <p>140225</p>

BECKER



For decades Becker quality steel furniture and equipment have been a fixed item in medical technology in clinics and hospitals.

Design, individual choice of colours as well as an extensive range of items and sizes allow you to make optimum use of your available space. The resulting concept is customized to your needs and logical through and through.

Metallwarenfabrik
Walter H. Becker GmbH
Anzenkirchener Strasse 4
84371 Triftern, Germany

Phone +49(0)8562 9605-0
Fax +49(0)8562 9605-64
info@becker-triftern.de
www.becker-triftern.de

Becker equipment:

- for the ward
- for the operating theatre
- for the central sterilisation unit
- for all functional rooms